



## NEBULIZER KIT

### DESCRIPTION

- Designed to be used with any standard Nebulizer kit.
- Super smooth kink resistance tubing ensures uniform flowrate.
- Smooth and feathered edge of face mask for patient comfort and reducing irritation points.
- Adjustable elastic strip and integrated nose clip for proper positioning of mask.
- The Nebulizer Kit just plugs into the standard fitting located at the bottom of the Mask.
- Masks are designed to deliver more medication directly to the patient due to the direct flow design.
- Proximal end of tube is connected with funnel shape connector for easy connection with oxygen source.
- Star lumen main tube to avoid accidental blockage.
- Tube Length : 210cm.
- Size : Adult, Paediatric & Neonatal.