



BIOHAZARD SAFETY BOX

The safety box is designed for safe disposal of medical sharps and needles which may otherwise result in injuries. The safety box is one of the most innovatively designed products by HMD healthcare and is a highly sought after medical consumable.

PRODUCT FEATURES:

- Special and Square shape
- Leak Proof and Puncture Resistant
- Volumes is available
- Double Layer corrugated paper board materials
- The medical waste sharp boxes are designed to permit safe disposal of sharps and medical waste. They can collect, store, transport and dispose of medical waste in hospital, clinic and family care
- Bio-degradable board material with Lifting Handle.
- Safe aperture to receive Sharp Waste.
- Possible to close and seal the aperture any time during use.
- Maximum filling line marking
- Stable, spill proof & damage resistant.
- Pictorial instructions to use.
- Available in 5 Ltrs., 10 Ltrs. and 20 Ltrs. Capacity
- As per WHO Standards: WHO/PQS/E10/SB01.1 S