



## GUEDEL AIRWAY

### INTENDED USE:

Guedel Airway is also known as an Oral Airway. It is an airway adjunct used to maintain or open the airway by stopping the tongue from covering the epiglottis.

### FEATURES:

- The Guedel Airway is manufactured from kink resistant and non-toxic polyethylene.
- The rounded edges are designed to minimize tissue damage.
- Bite block prevents biting of tongue and obstruction of airway.
- Consists of smooth airway path which aid in easy cleaning.
- Color coded for easy identification of sizes.
- Individually packed.
- Sterile, single use.

### SIZES

SIZES	COLOUR CODE	SIZES	COLOUR CODE
		2	Green
00	Blue	3	Yellow
0	Black	4	Red
1	White	5	Orange