



COMBINE DRESSING

Description:

The combine dressing is used to pack large bleeding wounds after initial bleeding control by direct pressure. They are wound dressings consisting of a highly absorbent thick layer of fleece enclosed in a soft and conformable non-woven fabric, folded over to form a flat tube. The non-woven covered pads have soft end seals to prevent pressure areas from forming.

Specification:

Usage/Application	Protection of Wound
Dressing Size	It is available in 10 X 20 cm and 10 X 10 cm pad and also roll of 10 or 20 cm wide X 3 mtr length
Material	Absorbent Cotton and Gauze
Type	Pad and Roll
Size	It is available in different sizes of Pad and Roll.
Color	White
Feature	Quick Healing