



## HIGH CONCENTRATION OXYGEN MASK

### Description:

- Super smooth kink resistance tubing ensures uniform flowrate.
- Smooth and feathered edge of face mask for patient comfort and reducing irritation points.
- High concentration Oxygen mask with tubing is intended to supply oxygen or other gases to an Individual for oxygen therapy, and the reservoir bag is able to avoid rebreathing.
- Proximal end of tube is connected with funnel shape connector for easy connection with oxygen source.
- Star lumen main tube to avoid accidental blockage.
- Tube Length: 210cm.
- Size: Adult & Paediatric.
- Available with DEHP free Material. Sterile / Disposable / Individually Packed.